



# Easy to use, Easy to handle

advanced features made simple

## quick tip

THE DVLA RECOMMENDS PATIENTS  
TEST BEFORE THEY DRIVE!



# TRUEresult



in partnership



# TRUEresult

easy to use, easy to handle  
advanced features made simple

## Advanced performance no-coding meter for easy, accurate testing at home.

- **No Coding** - eliminates a step that can lead to inaccurate results
- **Results As Fast As 4 Seconds** - means hardly any waiting time.
- **Tiny, 0.5 Microlitre Sample Size** - less blood means less painful testing
- **500 Test Memory** - includes time/date and 7-, 14- and 30-day averaging
- **Four Testing Reminder Alarms** - help promote testing compliance
- **Weekday Display** - for easy reference of test results by day
- **Ketone Test Reminder** - prompts testing for ketone levels as appropriate
- **Strip Release Button** - eliminates handling of strips after testing



Nipro Diagnostics (UK) Limited  
25 Barnes Wallis Road, Segensworth East  
Fareham, Hampshire, PO15 5TT  
**free meter upgrade ☎ 0800 08 588 08**

© 11/12 NIPRO DIAGNOSTICS LTD TR121005 11/12  
Reference 11 [www.dft.gov.uk/dvla/medical/ata glance.aspx](http://www.dft.gov.uk/dvla/medical/ata glance.aspx) - Last accessed Oct 12



in partnership

